

Protein disulfide-isomerase A3

Catalog No: #AP86384



Package Size: #AP86384-1 50ug #AP86384-2 100ug #AP86384-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Protein disulfide-isomerase A3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Endoplasmic reticulum resident protein 57,ER protein 57,Glucose-regulated thiol oxidoreductase 58 kDa protein,ERP57,GRP58
Accession No.	Uniprot ID:Q8JG64
Uniprot	Q8JG64
GeneID	373899;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PDIA3

Note: This product is for in vitro research use only