

Pyridoxal kinase

Catalog No: #AP86385



Package Size: #AP86385-1 50ug #AP86385-2 100ug #AP86385-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Pyridoxal kinase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Pyridoxine kinase
Accession No.	Uniprot ID:Q0II59
Uniprot	Q0II59
GeneID	514168;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: PDXK

Note: This product is for in vitro research use only