

P2Y purinoceptor 12

Catalog No: #AP86388



Package Size: #AP86388-1 50ug #AP86388-2 100ug #AP86388-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	P2Y purinoceptor 12
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	P2Y12,ADP-glucose receptor,ADPG-R,P2T(AC),P2Y(AC),P2Y(cyc),P2Y12 platelet ADP receptor,P2Y(ADP),SP1999,HORK3
Accession No.	Uniprot ID:Q9H244
Uniprot	Q9H244
GeneID	64805;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:P2RY12

Note: This product is for in vitro research use only