

# Putative phospholipase B-like 2

Catalog No: #AP86392

Package Size: #AP86392-1 50ug #AP86392-2 100ug #AP86392-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Putative phospholipase B-like 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	66.3 kDa protein,76 kDa protein,p76,LAMA-like protein 2,Lamina ancestor homolog 2,Phospholipase B domain-containing protein 2
Accession No.	Uniprot ID:Q3TCN2
Uniprot	Q3TCN2
GeneID	71772;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Plbd2

Note: This product is for in vitro research use only