

1-phosphatidylinositol 4,5-bisphosphate phosphodiesterase gamma-1

Catalog No: #AP86393

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP86393-1 50ug #AP86393-2 100ug #AP86393-3 1mg

Description

Product Name	1-phosphatidylinositol 4,5-bisphosphate phosphodiesterase gamma-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	PLC-148,Phosphoinositide phospholipase C-gamma-1,Phospholipase C-II,PLC-II,Phospholipase C-gamma-1,PLC-gamma-1,PLC1
Accession No.	Uniprot ID:P19174
Uniprot	P19174
GeneID	5335;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PLCG1

Note: This product is for in vitro research use only