

# 1-phosphatidylinositol 4,5-bisphosphate phosphodiesterase beta-1

Catalog No: #AP86395

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP86395-1 50ug #AP86395-2 100ug #AP86395-3 1mg

## Description

Product Name	1-phosphatidylinositol 4,5-bisphosphate phosphodiesterase beta-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	PLC-154, Phosphoinositide phospholipase C-beta-1, Phospholipase C-I, PLC-I, Phospholipase C-beta-1, PLC-beta-1, KIAA0581
Accession No.	Uniprot ID:Q9NQ66
Uniprot	Q9NQ66
GeneID	23236;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: PLCB1

Note: This product is for in vitro research use only