

## Pyridoxal phosphate phosphatase

Catalog No: #AP86396



Package Size: #AP86396-1 50ug #AP86396-2 100ug #AP86396-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Pyridoxal phosphate phosphatase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	PLP phosphatase,Chronophin,Reg I-binding protein 1,Cin,Plp,Plpp,Rbp1
Accession No.	Uniprot ID:Q8VD52
Uniprot	Q8VD52
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Pdxp

Note: This product is for in vitro research use only