

Phospholipase A2

Catalog No: #AP86401



Package Size: #AP86401-1 50ug #AP86401-2 100ug #AP86401-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Phospholipase A2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Group IB phospholipase A2, Phosphatidylcholine 2-acylhydrolase 1B, PLA2, PLA2A, PPLA2
Accession No.	Uniprot ID:P04054
Uniprot	P04054
GeneID	5319;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: PLA2G1B

Note: This product is for in vitro research use only