

Platelet-activating factor acetylhydrolase 2, cytoplasmic

Catalog No: #AP86402

Package Size: #AP86402-1 50ug #AP86402-2 100ug #AP86402-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Platelet-activating factor acetylhydrolase 2, cytoplasmic
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Serine-dependent phospholipase A2,SD-PLA2
Accession No.	Uniprot ID:P83006
Uniprot	P83006
GeneID	313611;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Pafah2

Note: This product is for in vitro research use only