

Plasminogen activator inhibitor 1

Catalog No: #AP86403



Package Size: #AP86403-1 50ug #AP86403-2 100ug #AP86403-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Plasminogen activator inhibitor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	PAI,Endothelial plasminogen activator inhibitor,Serpin E1,Pai1,Planh1
Accession No.	Uniprot ID:P20961
Uniprot	P20961
GeneID	24617;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Serpine1

Note: This product is for in vitro research use only