

Metaxin-2

Catalog No: #AP86407



Package Size: #AP86407-1 50ug #AP86407-2 100ug #AP86407-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Metaxin-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Mitochondrial outer membrane import complex protein 2
Accession No.	Uniprot ID:Q2L969
Uniprot	Q2L969
GeneID	733591;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: MTX2

Note: This product is for in vitro research use only