

Myogenic factor 6

Catalog No: #AP86408



Package Size: #AP86408-1 50ug #AP86408-2 100ug #AP86408-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Myogenic factor 6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Myf-6, Muscle-specific regulatory factor 4, MRF4
Accession No.	Uniprot ID:Q7YS80
Uniprot	Q7YS80
GeneID	281336;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: MYF6

Note: This product is for in vitro research use only