

CYB561D1 Antibody

Catalog No: #44313

Package Size: #44313-1 50ul #44313-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	CYB561D1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	CYB561D1 Antibody detects endogenous levels of total CYB561D1
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human CYB561D1
Target Name	CYB561D1
Other Names	CYB561D1, RP5-831G13.3
Accession No.	Swiss-Prot#: Q8N8Q1NCBI Gene ID:284613
Uniprot	Q8N8Q1
GeneID	284613;
Calculated MW	23kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C


Application Details

Western blotting: 1:1000

Images



Western blot analysis of lysates from HeLa cells and COLO205 cells, using Cytochrome b561 D1 Antibody. The lane on the right is blocked with the synthesized peptide.



Immunofluorescence analysis of MCF-7 cells, using Cytochrome b561 D1 Antibody. The picture on the right is blocked with the synthesized peptide.

Note: This product is for in vitro research use only