

Myosin light chain kinase, smooth muscle

Catalog No: #AP86410



Package Size: #AP86410-1 50ug #AP86410-2 100ug #AP86410-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Myosin light chain kinase, smooth muscle
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	MLCK,Telokin
Accession No.	Uniprot ID:Q28824
Uniprot	Q28824
GeneID	338037;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MYLK

Note: This product is for in vitro research use only