

Serine/threonine-protein kinase PLK1

Catalog No: #AP86413



Package Size: #AP86413-1 50ug #AP86413-2 100ug #AP86413-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Serine/threonine-protein kinase PLK1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Polo-like kinase 1,PLK-1,Plk
Accession No.	Uniprot ID:Q62673
Uniprot	Q62673
GeneID	25515;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Plk1

Note: This product is for in vitro research use only