

Dolichyl-diphosphooligosaccharide--protein glycosyltransferase 48 kDa subunit

Catalog No: #AP86414

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP86414-1 50ug #AP86414-2 100ug #AP86414-3 1mg

Description

| | |
|--------------------|---|
| Product Name | Dolichyl-diphosphooligosaccharide--protein glycosyltransferase 48 kDa subunit |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Rt |
| Other Names | DDOST 48 kDa subunit |
| Accession No. | Uniprot ID:Q641Y0 |
| Uniprot | Q641Y0 |
| GeneID | 313648; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ddost

Note: This product is for in vitro research use only