

Beta-microseminoprotein

Catalog No: #AP86415



Package Size: #AP86415-1 50ug #AP86415-2 100ug #AP86415-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Beta-microseminoprotein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Immunoglobulin-binding factor,IGBF,PN44,Prostate secreted seminal plasma protein,Prostate secretory protein of 94 amino acids,PSP-94,Seminal plasma beta-inhibin,PRSP
Accession No.	Uniprot ID:P08118
Uniprot	P08118
GeneID	4477;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MSMB

Note: This product is for in vitro research use only