

Pro-neuregulin-1, membrane-bound isoform

Catalog No: #AP86416



Package Size: #AP86416-1 50ug #AP86416-2 100ug #AP86416-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Pro-neuregulin-1, membrane-bound isoform
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Pro-NRG1,ARIA
Accession No.	Uniprot ID:Q05199
Uniprot	Q05199
GeneID	373906;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:NRG1

Note: This product is for in vitro research use only