

Phospholipase A2, membrane associated

Catalog No: #AP86418



Package Size: #AP86418-1 50ug #AP86418-2 100ug #AP86418-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Phospholipase A2, membrane associated
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Enhancing factor,EF,GIIC sPLA2,Group IIA phospholipase A2,Phosphatidylcholine 2-acylhydrolase 2A
Accession No.	Uniprot ID:P31482
Uniprot	P31482
GeneID	18780;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Pla2g2a

Note: This product is for in vitro research use only