

Pituitary adenylate cyclase-activating polypeptide

Catalog No: #AP86419



Package Size: #AP86419-1 50ug #AP86419-2 100ug #AP86419-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Pituitary adenylate cyclase-activating polypeptide
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	PACAP
Accession No.	Uniprot ID:P13589
Uniprot	P13589
GeneID	24166;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Adcyap1

Note: This product is for in vitro research use only