

## CYP4F2 Antibody

Catalog No: #44314

Package Size: #44314-1 50ul #44314-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	CYP4F2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB;IHC;IF
Species Reactivity	Hu
Specificity	CYP4F2 Antibody detects endogenous levels of total CYP4F2
Immunogen Type	Peptide
Immunogen Description	The antiserum was produced against synthesized peptide derived from human CYP4F2.
Target Name	CYP4F2
Other Names	CYP4F2, CYP1VF2, Cytochrome P450 4F2, Cytochrome P450-LTB-omega, CPF2
Accession No.	Swiss-Prot#: P78329NCBI Gene ID:8529
Uniprot	P78329
GeneID	8529;
Calculated MW	60kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C


## Application Details

WB 1:500 - 1:2000. IHC 1:100 - 1:300. IF 1:50-200

## Images



Western Blot analysis of HeLa cells using CYP4F2 Polyclonal Antibody



Immunohistochemistry analysis of paraffin-embedded human liver carcinoma, using Cytochrome P450 4F2 Antibody. The picture on the right is blocked with the synthesized peptide.

---

Note: This product is for in vitro research use only