

Sodium/calcium exchanger 1

Catalog No: #AP86424



Package Size: #AP86424-1 50ug #AP86424-2 100ug #AP86424-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Sodium/calcium exchanger 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Na(+)/Ca(2+)-exchange protein 1,Solute carrier family 8 member 1,Ncx1
Accession No.	Uniprot ID:Q01728
Uniprot	Q01728
GeneID	29715;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Slc8a1

Note: This product is for in vitro research use only