

N-acyl-phosphatidylethanolamine-hydrolyzing phospholipase D

Catalog No: #AP86425

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP86425-1 50ug #AP86425-2 100ug #AP86425-3 1mg

Description

Product Name	N-acyl-phosphatidylethanolamine-hydrolyzing phospholipase D
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	N-acyl phosphatidylethanolamine phospholipase D,Mblcd1
Accession No.	Uniprot ID:Q8BH82
Uniprot	Q8BH82
GeneID	242864;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Napepld

Note: This product is for in vitro research use only