

Oxidized low-density lipoprotein receptor 1

Catalog No: #AP86426



Package Size: #AP86426-1 50ug #AP86426-2 100ug #AP86426-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Oxidized low-density lipoprotein receptor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Ox-LDL receptor 1,C-type lectin domain family 8 member A,Lectin-like oxidized LDL receptor 1,LOX-1,Lectin-type oxidized LDL receptor 1,CLEC8A,LOX1
Accession No.	Uniprot ID:P78380
Uniprot	P78380
GeneID	4973;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:OLR1

Note: This product is for in vitro research use only