

# Myosin-1

Catalog No: #AP86428



Package Size: #AP86428-1 50ug #AP86428-2 100ug #AP86428-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Myosin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Myosin heavy chain 1,Myosin heavy chain 2x,MyHC-2x,Myosin heavy chain, skeletal muscle, adult 1
Accession No.	Uniprot ID:Q076A6
Uniprot	Q076A6
GeneID	100135061;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:MYH1

Note: This product is for in vitro research use only