

Ornithine aminotransferase, mitochondrial

Catalog No: #AP86430



Package Size: #AP86430-1 50ug #AP86430-2 100ug #AP86430-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ornithine aminotransferase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Ornithine delta-aminotransferase, Ornithine--oxo-acid aminotransferase
Accession No.	Uniprot ID:P04181
Uniprot	P04181
GeneID	4942;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: OAT

Note: This product is for in vitro research use only