

Beta-type platelet-derived growth factor receptor

Catalog No: #AP86433



Package Size: #AP86433-1 50ug #AP86433-2 100ug #AP86433-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Beta-type platelet-derived growth factor receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	PDGF-R-beta,CD140 antigen-like family member B,Beta-type platelet-derived growth factor receptor,CD140 antigen-like family member B,Platelet-derived growth factor receptor 1,PDGFR-1,CD140b,PDGFR,PDGFR1
Accession No.	Uniprot ID:Q6QNF3
Uniprot	Q6QNF3
GeneID	442985;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PDGFRB

Note: This product is for in vitro research use only