

Prolyl 4-hydroxylase subunit alpha-2

Catalog No: #AP86436



Package Size: #AP86436-1 50ug #AP86436-2 100ug #AP86436-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Prolyl 4-hydroxylase subunit alpha-2 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Chicken |
| Other Names | 4-PH alpha-2,Procollagen-proline,2-oxoglutarate-4-dioxygenase subunit alpha-2,RCJMB04_5117 |
| Accession No. | Uniprot ID:Q5ZLK5 |
| Uniprot | Q5ZLK5 |
| GeneID | 416326; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:P4HA2

Note: This product is for in vitro research use only