

Phospholipase A2

Catalog No: #AP86437



Package Size: #AP86437-1 50ug #AP86437-2 100ug #AP86437-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Phospholipase A2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Group IB phospholipase A2, PLA2-Ib, Phosphatidylcholine 2-acylhydrolase 1B, Pla2
Accession No.	Uniprot ID:Q9Z0Y2
Uniprot	Q9Z0Y2
GeneID	18778;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Pla2g1b

Note: This product is for in vitro research use only