

Osteopontin

Catalog No: #AP86438



Package Size: #AP86438-1 50ug #AP86438-2 100ug #AP86438-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Osteopontin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Bone sialoprotein 1, Secreted phosphoprotein 1, SPP-1
Accession No.	Uniprot ID:P23498
Uniprot	P23498
GeneID	395210;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: SPP1

Note: This product is for in vitro research use only