

Retrotransposon-derived protein PEG10

Catalog No: #AP86445



Package Size: #AP86445-1 50ug #AP86445-2 100ug #AP86445-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Retrotransposon-derived protein PEG10
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Embryonal carcinoma differentiation-regulated protein, Mammalian retrotransposon-derived protein 2, Myelin expression factor 3-like protein 1, MEF3-like protein 1, Paternally expressed gene 10 protein, Retrotransposon gag domain-containing protein 3, Retrotransposon-derived gag-like polyprotein, Ty3/Gypsy-like protein, EDR, KIAA1051, MAR2, MART2, MEF3L1, RGAG3
Accession No.	Uniprot ID: Q86TG7
Uniprot	Q86TG7
GeneID	23089;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: PEG10

Note: This product is for in vitro research use only