

High affinity nerve growth factor receptor

Catalog No: #AP86453



Package Size: #AP86453-1 50ug #AP86453-2 100ug #AP86453-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	High affinity nerve growth factor receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Neurotrophic tyrosine kinase receptor type 1
Accession No.	Uniprot ID:Q3UFB7
Uniprot	Q3UFB7
GeneID	18211;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Ntrk1

Note: This product is for in vitro research use only