

Proteasome subunit beta type-9

Catalog No: #AP86455



Package Size: #AP86455-1 50ug #AP86455-2 100ug #AP86455-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Proteasome subunit beta type-9
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Proteasome subunit beta-1i
Accession No.	Uniprot ID:Q3SZC2
Uniprot	Q3SZC2
GeneID	510593;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: PSMB9

Note: This product is for in vitro research use only