

## PC4 and SFRS1-interacting protein

Catalog No: #AP86456



Package Size: #AP86456-1 50ug #AP86456-2 100ug #AP86456-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	PC4 and SFRS1-interacting protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CLL-associated antigen KW-7,Dense fine speckles 70 kDa protein,DFS 70,Lens epithelium-derived growth factor,Transcriptional coactivator p75/p52,DFS70,LEDGF,PSIP2
Accession No.	Uniprot ID:O75475
Uniprot	O75475
GeneID	11168;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:PSIP1

Note: This product is for in vitro research use only