

Proteinase-activated receptor 4

Catalog No: #AP86460



Package Size: #AP86460-1 50ug #AP86460-2 100ug #AP86460-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Proteinase-activated receptor 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	PAR-4,Coagulation factor II receptor-like 3,Thrombin receptor-like 3,PAR4
Accession No.	Uniprot ID:Q96RI0
Uniprot	Q96RI0
GeneID	9002;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:F2RL3

Note: This product is for in vitro research use only