

# Phosphatidylinositol 4,5-bisphosphate 3-kinase catalytic subunit delta isoform

Catalog No: #AP86462

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP86462-1 50ug #AP86462-2 100ug #AP86462-3 1mg

## Description

Product Name	Phosphatidylinositol 4,5-bisphosphate 3-kinase catalytic subunit delta isoform
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	PI3-kinase subunit delta,Phosphatidylinositol 4,5-bisphosphate 3-kinase 110 kDa catalytic subunit delta,PtdIns-3-kinase subunit p110-delta
Accession No.	Uniprot ID:O00329
Uniprot	O00329
GeneID	5293;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:PIK3CD

Note: This product is for in vitro research use only