

Decorin

Catalog No: #AP86464



Package Size: #AP86464-1 50ug #AP86464-2 100ug #AP86464-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Decorin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Bone proteoglycan II,PG-S2
Accession No.	Uniprot ID:Q9XSD9
Uniprot	Q9XSD9
GeneID	396957;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:DCN

Note: This product is for in vitro research use only