

# 2-oxoglutarate dehydrogenase, mitochondrial

Catalog No: #AP86466



Package Size: #AP86466-1 50ug #AP86466-2 100ug #AP86466-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	2-oxoglutarate dehydrogenase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	2-oxoglutarate dehydrogenase complex component E1,OGDC-E1,Alpha-ketoglutarate dehydrogenase,Kiaa4192
Accession No.	Uniprot ID:Q60597
Uniprot	Q60597
GeneID	18293;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Ogdh

Note: This product is for in vitro research use only