

# Oligodendrocyte transcription factor 1

Catalog No: #AP86467



Package Size: #AP86467-1 50ug #AP86467-2 100ug #AP86467-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Oligodendrocyte transcription factor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Oligo1,Olg-1 bHLH protein
Accession No.	Uniprot ID:Q9WUQ3
Uniprot	Q9WUQ3
GeneID	60394;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Olig1

Note: This product is for in vitro research use only