

Neurotrimin

Catalog No: #AP86468



Package Size: #AP86468-1 50ug #AP86468-2 100ug #AP86468-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Neurotrimin |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Rt |
| Other Names | GP65,Hnt,Nt |
| Accession No. | Uniprot ID:Q62718 |
| Uniprot | Q62718 |
| GeneID | 50864; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ntm

Note: This product is for in vitro research use only