

Platelet endothelial cell adhesion molecule

Catalog No: #AP86472



Package Size: #AP86472-1 50ug #AP86472-2 100ug #AP86472-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Platelet endothelial cell adhesion molecule
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	PECAM-1,CD31,Pecam,Pecam-1
Accession No.	Uniprot ID:Q08481
Uniprot	Q08481
GeneID	18613;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Pecam1

Note: This product is for in vitro research use only