

## Polycomb group RING finger protein 6

Catalog No: #AP86476

Package Size: #AP86476-1 50ug #AP86476-2 100ug #AP86476-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Polycomb group RING finger protein 6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	RING finger protein 134,Rnf134
Accession No.	Uniprot ID:Q5XI70
Uniprot	Q5XI70
GeneID	309457;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Pcgf6

Note: This product is for in vitro research use only