

Olfactory receptor 51E2

Catalog No: #AP86477



Package Size: #AP86477-1 50ug #AP86477-2 100ug #AP86477-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Olfactory receptor 51E2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	G-protein coupled receptor RA1c,Olfactory receptor 59,Olr59,Psg
Accession No.	Uniprot ID:O88628
Uniprot	O88628
GeneID	170816;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Or51e2

Note: This product is for in vitro research use only