

## Protein phosphatase 1 regulatory subunit 14B

Catalog No: #AP86481



Package Size: #AP86481-1 50ug #AP86481-2 100ug #AP86481-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Protein phosphatase 1 regulatory subunit 14B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Phosphatase holoenzyme inhibitor 1, PHI-1, Ubiquitous PKC-potentiated PP1 inhibitor, Phi1
Accession No.	Uniprot ID:Q8K3F3
Uniprot	Q8K3F3
GeneID	259225;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Ppp1r14b

Note: This product is for in vitro research use only