

Peptidyl-prolyl cis-trans isomerase-like 1

Catalog No: #AP86482



Package Size: #AP86482-1 50ug #AP86482-2 100ug #AP86482-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Peptidyl-prolyl cis-trans isomerase-like 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	PPlase, Rotamase PPIL1, CYPL1, CGI-124, UNQ2425/PRO4984
Accession No.	Uniprot ID:Q9Y3C6
Uniprot	Q9Y3C6
GeneID	51645;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: PPIL1

Note: This product is for in vitro research use only