

Nucleoredoxin-like protein 1

Catalog No: #AP86483



Package Size: #AP86483-1 50ug #AP86483-2 100ug #AP86483-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Nucleoredoxin-like protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Thioredoxin-like protein 6, TXNL6
Accession No.	Uniprot ID:Q96CM4
Uniprot	Q96CM4
GeneID	115861;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: NXNL1

Note: This product is for in vitro research use only