

Pyridoxal kinase

Catalog No: #AP86485



Package Size: #AP86485-1 50ug #AP86485-2 100ug #AP86485-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Pyridoxal kinase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Pyridoxine kinase,PKH
Accession No.	Uniprot ID:O46560
Uniprot	O46560
GeneID	396983;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PDXK

Note: This product is for in vitro research use only