

Myeloperoxidase

Catalog No: #AP86486



Package Size: #AP86486-1 50ug #AP86486-2 100ug #AP86486-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Myeloperoxidase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	MPO
Accession No.	Uniprot ID:P11247
Uniprot	P11247
GeneID	17523;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Mpo

Note: This product is for in vitro research use only