

NADPH oxidase 4

Catalog No: #AP86487



Package Size: #AP86487-1 50ug #AP86487-2 100ug #AP86487-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NADPH oxidase 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Kidney oxidase-1,KOX-1,Kidney superoxide-producing NADPH oxidase,Renal NAD(P)H-oxidase,Superoxide-generating NADPH oxidase 4,Renox
Accession No.	Uniprot ID:Q9JHI8
Uniprot	Q9JHI8
GeneID	50490;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Nox4

Note: This product is for in vitro research use only