

POU domain, class 3, transcription factor 1

Catalog No: #AP86488



Package Size: #AP86488-1 50ug #AP86488-2 100ug #AP86488-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	POU domain, class 3, transcription factor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Octamer-binding protein 6, Oct-6, Octamer-binding transcription factor 6, OTF-6, POU domain transcription factor SCIP, OCT6, OTF6
Accession No.	Uniprot ID: Q03052
Uniprot	Q03052
GeneID	5453;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: POU3F1

Note: This product is for in vitro research use only